

Adjustable Height Rollator, 6" Casters

HCPCS: E0143+E0156



FIGURE A



FIGURE B



FIGURE C




FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Universal height adjustment on frame allows seat height to be adjusted from 18" to 24" in 1" increments
- Comes with a comfortable, durable padded seat
- Zippered pouch under seat for added privacy and security of personal items
- Removable, hinged, padded backrest can be folded up or down
- 6" casters with soft grip tires ideal for indoor and outdoor use and make the ride smooth and stable
- Ergonomic handles are easy-to-grip and are adjustable in height with self-threading knob

PRODUCTS

Description

RTL10261BL	ADJ Height Rollator,6" Wheels,Blue,1EA
RTL10261RD	ADJ Height Rollator,6" Wheels,Red,1EA

SPECIFICATIONS

Description

Carton Shipping Weight	16 lbs.
Casters	6"
Handle Height	29.5" - 38"
Inside Hand Grip Width	17"
Number of Wheels	4
Seat Depth	12"
Seat To Floor Height	18" - 22"
Seat Width	13"
Unit Dimensions	24" (W) x 27" (L)
Warranty	Limited Lifetime
Weight	14 lbs.
Weight Capacity	300 lbs

PARTS

Description

10216-01	Brake Assembly, 10216, 1/ea, EA
9502M773412	Front Wheel, 7734 W/ 2M Serial #, 1/ea, EA
9505W1026102	Height Adj Knob,10261,1pair, PR
9505W1026105	Folding Hinge Assembly,10261,1ea, EA
9505W1026106	Right Rear Leg w/Axle Sleeve,10261,1ea, EA
9505W1026107	Left Rear Leg w/Axle Sleeve,10261,1ea, EA
9505W1026108	Right Wheel Brake,10261,1ea, EA
9505W1026109	Left Wheel Brake,10261,1ea, EA
9505W1026111	Rear Axle Bolt Nut Washers,10261,1ea, EA
9505W1026110	Rear Wheel, 10261, 1/ea, EA

PARTS	Description
9505W1026113	Front Axle Nut Washers,10261,1ea, EA
9505W1026114	Front Fork,10261,1ea, EA
9505W1026104	Seat Hinge Brackets,10261,1pair, PR
9505W1026115	Front Leg,Left,10261,1ea, EA
9505W1026123	Leg Pin for 10261,4/bag, BAG
9505W1026122	Frame Plug for 10261,4/bag, BAG
9505W1026121	BackrestPivotBracket10261,1pr, PR
9505W1026119	Handle Bar,Right,10261,1ea, EA
9505W1026118	Handle Bar,Left,10261,1ea, EA
9501S72603	Backrest Pivot Brackets for R726,1 pr, PR