

Deluxe Bariatric Shower Chair with Cross-Frame Brace

HCPCS: E0245



FIGURE A

- Frame has width/depth cross braces attached with aircraft-type rivets
- Large seat with 16.5" depth
- Seat height adjusts in 1" increments
- Durable, blow-molded plastic bench
- Tool-free removable back
- Aluminum frame is lightweight, durable and corrosion proof
- Drainage holes reduce slipping
- Weight Capacity 600 lbs
- Made using recycled materials

PRODUCTS

Description

12029-2	BathBenchOverszd.HD.w/Bk,2/cs
---------	-------------------------------

SPECIFICATIONS

Description

Carton Shipping Weight	10 20 lbs.
Outside Legs	21.25" (W) x 18.5" (D)
Seat Dimensions	16.5" (W) x 16.5" (D)
Seat Height	16.5" - 21.5"
Warranty	Limited Lifetime
Weight	9 lbs.
Weight Capacity	600 lbs

ACCESSORIES

Description

12028T	Tips, BathBenchOverSized,12029, HD, 4/bx, BX
--------	--