

Gel-U-Seat Lite General Use 2" Gel/Foam Wheelchair Cushion

HCPCS: E2601



FIGURE A



FIGURE B

- Offering great versatility, the Gel-U-Seat Lite provides comfort to almost all wheelchair users
- Unique dual-chambered gel bladder keeps the pelvis in contact with pressure redistributing cooling gel
- The cover is fluid resistant low-shear nylon zippered and removable with a durable vinyl bottom
- Secures easily to wheelchair with 12" tie straps
- Meets flammability standard TB117

PRODUCTS

Description

8040-1	2" Gel Cushion 16x16x2, 1 ea,4/CS
8040-2	2" Gel Cushion 18x16x2, 1 ea,4/CS
8040-3	2" Gel Cushion 20x16x2, 1 ea,4/CS
8040-4	2" Gel Cushion 20x18x2, 1 ea,4/CS
8040-5	2" Gel Cushion 18x18x2", 1 ea,4/CS
8040-7	2" Gel Cushion 22x18x2", 1 ea/4CS
8040-9	2" Gel Cushion 24x18x2", 1 ea,4/CS

SPECIFICATIONS

Description

Cover Material	Fluid-Resistant Stretch Top Waterproof Vinyl Bottom
Warranty	1 Year Limited
Weight Capacity	Up to 20" 275 lbs, 22" and up 500 lbs