


# AL100 / AL100HD

True Blue Outside Lift You Can Count On

**\$3400.00 INSTALLED**



920 TATE BLVD SE SUITE 104  
HICKORY, NC 28602  
828-855-9099  
www.thewheelchairplace.com



## A DURABLE POWERHOUSE

The AL100 / AL100HD are durable lifters with weight capacities of 350 lb / 400 lb. The AL100HD is the heavy lifter and is the only 400 lb weight capacity on the market. These universal scooter lifts are the most popular on the market. Why? Because of its durable aluminum platform, engineering, and effortless operation, this lift endures the test of time and weather. Its unique automatic hold-down foot secures the scooter. No scooter modifications are required for this lift.

We design our lifts with you in mind. Offering more customized platform sizes that are 15 to 25% lighter than others, which means you save money on fuel.

## EASY CARGO ACCESS

Our lightweight swing-away is our latest accessory for your outside auto lift. This innovative design uses a new I-beam construction, resulting in a 13 lb reduction that provides you with greater handling and steering.

# YOU CAN COUNT ON IT


AL100/AL100HD



Toggle switch to rise and lower lift with ease.



Automatic hold-down foot secures your scooter and power chair for you. No bending needed.



Harmer's patriotic Stars -N- Stripes platforms are engineered for maximum strength.



*I've had my AL100 for years. It's true BLUE, easy to use, and stands up against the harshest of climates.*

## 83

Pound Install Weight

## 350

Pound Capacity

## 84

Pound Install Weight

## 400

Pound Capacity

### AL100 Product Specifications

Accommodates	Scoters
Platform	26" x 39" Aluminum Deck; Adjustable Cradle Extends up to 46"
Lifting Capacity	350 lb
Installed Weight	83 lb
Securement	Automatic hold-down Foot and Rear Wheel Cradle
Hitch	II or greater
Options	Battery Pack, Swing-Away, Scooter Covers
Warranty	3-Year Transferable

### AL100HD Product Specifications

Accommodates	Scoters
Platform	26" x 39" Aluminum Deck; Adjustable Cradle Extends up to 46"
Lifting Capacity	400 lb
Installed Weight	84 lb
Securement	Automatic hold-down Foot and Rear wheel Cradle
Hitch	III or greater
Options	Battery Pack, Swing-Away, Scooter Covers
Warranty	3-Year Transferable

AL100 / 02110218



Harmar.com | 800.833.0478

