



920 TATE BLVD SE SUITE 104
HICKORY, NC 28602
828-855-9099
www.thewheelchairplace.com



Innovative, Lightweight Wheelchair Lift

Lightweight, non-powered solution for manual folding wheelchairs. Innovative tilting joint helps the wheelchair roll on/off easily; adjustable hold-down arm safely secures it. Compact design folds down for tail gate / lift gate access, folds up when not needed.

Features & Benefits

- Tilting joint for roll-on ease
- Adjustable hold-down arm
- Folds to access tail gate / lift gate; folds up when not in use
- No chair modifications required
- Ball mount standard; upgrades available
- No hitch post needed with ball mount

Specifications	AL003 Tilt-N-Tote
Accommodates	Manual folding wheelchairs
Power	Manual
Lifting Capacity	100 lbs. (45 kg)
Installed Weight	25 lbs. (11.3 kg)
Securement	Manual hold-down arm
Hitch	Ball mount or greater
Standard Hitch Height	5" – 19"
Options	Manual wheelchair covers
Warranty	3-year transferable



Adjustable hold-down arm safely secures wheelchair for transport.

