INVACARE® 9000 TOPAZ™ WHEELCHAIR





The Invacare® 9000 Topaz™ Wheelchair is designed specifically to suit the unique needs of bariatric patients. The Topaz wheelchair combines the strength and durability of a carbon steel frame with the comfort of heavy-duty front riggings and arms for proper positioning. With a variety of seat widths and depths, frame reinforcements that include gussetted sides, caster journals and dual heavy-duty crossbraces, this chair provides not only ease in maneuverability but proven stability.

Features

Available in vinyl or nylon upholstery

Durable low-maintenance, carbon steel frame is long-lasting

Reinforced frame includes gussetted sides, caster journals and dual heavy-duty crossbraces

Multiple frame color variations available for customizing

Heavy-duty front riggings provide proper positioning for feet with footrest and elevating legrest

Heavy-duty arms designed to accommodate weight capacity and provide proper positioning

Dual axle positions allow for quick and simple seat-to-floor height variations, 17.5" and 19.5"

Heavy-duty package option for 1,000 lb. weight capacity

Heavy-duty triple inner liners keep seat and back from stretching

HCPCS Code: K0007













- 1. Dual axle positions
- 2. Heavy-duty front riggings
- 3. Heavy-duty arms
- 4. Bariatric sizes
- 5. Seating available
- 6. Builder options



Invacare Corporation

www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668 - 5324

©2012

Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 11-230 111247















Model # Overall Overall Overall Depth Folded Seat Seat Seat-to-Width Height w/o riggings Width Width Depth Floor Height (2" increments) (2" increments) 9TPZ 39" 32" 20"-30" 15.5"-19.5" 30" 16" 18", 20"

















Wheel Back Caster Weight Shipping Carton Limited Height Height Height Capacity Weight Dimensions Warranty 9TPZ 700 lb., 16"-18" 6", 8" 20", 24" 97 lb. 35"x15"x 46.5" 3 Years 1,000 lb



All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.