

Med-Aire Alternating Pressure Pump and Pad System

HCPCS: E0181



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Quiet pump alternately inflates and deflates 130 individual bubble air cells
- 4 LPM (liters per minute) pump produces consistent air flow and pressure over a 5 minute cycle time and is available fixed or variable (Fig.A)
- Variable pressure setting on pump allows for individualized comfort increasing patient compliance
- Durable, heavy gauge, premium vinyl bubble pad provides superior therapeutic treatment and is available with or without end flaps
- Pump has built-in brackets for mounting easily to bed

PRODUCTS

Description

14001E	APP Sys,Dlx,Variable Pressure,1/ea
14001EFD	Dlx APP w/Pad w/End Flaps,1/CV, 6/CT
14002E	APP Sys,Standard,Fixed Pressure,1/ea,RTL

SPECIFICATIONS

Description

Weight Capacity	300 lbs
-----------------	---------

ACCESSORIES

Description

14001T	Tubing for 14001/14002,1/ea
14003-EF	APP Pad, w/End Flaps,1/ea, RTL
14003	APP Pad Only, 1/ea, RTL
14005E	Pump only-14001E, BA9600-P, 1ea, EA
14006E	Pump only for 14002E,1/ea