

## 22" Bariatric Aluminum Transport Chair

HCPCS: E1039



FIGURE A



FIGURE B

- Weighs only 33 lbs
- Supports individuals up to 450 lbs
- Attachable push bar with padded foam grip (Figure A)
- 12" rear flat-free wheels
- Comes with seat belt
- Comes with swing-away footrests and anti-tippers
- Back folds down for convenient storage (Figure B)
- Padded armrests provide additional comfort
- Heavy-duty, nylon reinforced upholstery
- Comes with hand wheel lock lever and companion activated wheel locks

### PRODUCTS

### Description

ATC22-R	Aluminum HD 22" Red Transport Chair,1/ea
---------	--

### SPECIFICATIONS

### Description

Back Height	16"
Carton Shipping Weight	39 lbs
Overall Height	37.5"
Overall Length w/ Riggings	41"
Overall Width	30" (Open)   12" (Closed)
Seat Dimensions	22" (W) x 18" (D)
Seat To Floor Height	21"
Warranty	Limited Lifetime
Weight	33 lbs
Weight Capacity	450 lbs