

Premium Guard Gel Overlay

HCPCS: E0185



FIGURE A

- Best in class gel overlay optimizes pressure redistribution and is now available in up to 500 lbs weight capacity
- 15 gel bladders cover entire patient contact area for the prevention and treatment of pressure ulcers
- Full length "head-to-toe" gel bladders minimize "bottoming out" and "pitting"
- New convoluted foam topper vastly improves patient comfort and air circulation
- 4 elastic corner straps secure overlay to mattress
- Urethane coated nylon cover is water resistant and vapor permeable
- Vinyl base is stain resistant and waterproof
- Bladders contain an aqueous non-toxic gel and polystyrene beads to minimize gel migration

PRODUCTS

Description

14893	PremGuardGel/FoamOverlay34x76x3.5in,1/cs
14901	Prem.GuardGel/FoamOverlay42x76x3.5",1/cs

SPECIFICATIONS

Description

Weight Capacity	300-500 lbs
-----------------	-------------

ACCESSORIES

Description

14893C	Prem.Guard Gel,3"Foam Overlay,Cover,1/ea
14901C	Cover for 14901, 1/ea