

Med-Aire Plus 8" Alternating Pressure and Low Air Loss Mattress System

HCPCS: E0277



FIGURE A




FIGURE B

- Active combination therapy system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and microclimate control
- Designed to prevent, treat and heal pressure ulcers in the home or long term care setting
- Digital control unit offers 4 cycle times and adjustable comfort weight settings for optimal immersion and patient satisfaction
- 3 diagnostic audio and visual alarms alert for issues with air pressure, the system, or power
- System includes 20 individual air cells to maximize immersion and envelopment with pillow feature for head stability & comfort
- Advanced Clinical and Care Giver Features: Auto-Firm, Seat Inflate, microprocessing feedback, lock-out and more!
- Cell-on-cell designed air cells prevent "bottoming out" and provide up to 24 hours of power outage protection
- 8" deep air cells are easily removed with quick disconnects for convenient servicing
- Removable 4 way stretch PU cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- CPR valve allows for rapid deflation
- Tubing quickly and easily disconnects with EZ Lock/Quick Release design

PRODUCTS

Description

14029	MedAire8 LowAirLoss.w/Pump system,1/cs
14029-84	Med Aire 8" Low Air Loss APM w/pump 84

SPECIFICATIONS

Description

Weight Capacity	450 lbs
-----------------	---------

ACCESSORIES

Description

14427	Pump, Electric, for 14425,1/bx
-------	--------------------------------

PARTS

Description

14029XP	Digital 3-Outlet Pump 14029/14029DP,1/ea
14029-CELL	Bladder Only,14029,1/ea"
14029B	Carry bag for 14029,1/ea
14029MN	Mattress Only,14029,3Tube System,1/ea
14013	CPRPlugForLowAirLossMatt,3Port,1ea
14029-CELL-1	Cell Only wLazer Holes,14029,1/ea
DS16C01207	Mattress Hose 3 port, 1/ea
14029-FFTCPLUG	Plug For Connector,14029 Cell,Female,1e
14029C	Cover for 14029,1/ea
15005PC	Power Cord,1/ea

PARTS	Description
14001-L	L-Shape Connector, LAL,1/ea
14001-T	T-Connector for LAL, 1/ea
14029-FFTC	Connector for 14029 Cell, Female, 1/ea
14028V	OneWayValveForLowAirLossMatt,1/ea
14029-FMTC	Connector for 14029 Cell, Male, 1/ea
14028T	Tubing,7", LAL, 1/ea
14030XP-S	Pump Socket 3 Hose Connection,1/ea
14007-3PORT	Tubing Adapter 3-Port,14029DP,1/ea
14030HK	Hanging Hooks,14030XP,1pair
14029FN	Filter14029XP,1ea