



Experience the beauty of lightweight travel

with the Jazzy® Passport! This compact power chair folds in just a few simple steps for effortless transport. With convenient storage the Jazzy® Passport makes traveling a breeze.

Folded Size 31" x 16" x 23.5"









Weight Capacity 250 lbs. Maximum Speed¹ Up to 3.6 mph Ground Clearance³ 3" at motor Turning Radius³ 24.25"	
Ground Clearance ³ 3" at motor	
Turning Radius³ 24.25"	
Overall Length ³ 40"	
Overall Width³ 23.5"	
Seat-to-Ground Height 19"	
Seat-to-Foot Platform Height 14.5" (no adjustment)	
Casters 8" solid	
Drive Tires 12" Pneumatic drive	
Suspension Limited	
Range Per Charge ^{1,2} Up to 9.6 miles (standard) Up to 5.7 miles (optional)	
Total Weight Without Batteries ^{3,5} 54 lbs. (unit does not disassemble)	
Heaviest Piece When Disassembled ³ N/A (unit does not disassemble)	
Standard Seating Material: Gray, Foam Dimensions: 17.5" x 16.5"	
Drivetrain 2x 180W brushed motors	
Dual Braking System Regenerative and electromechanical	
Standard Electronics 50 amp, PG VSI	
Battery Requirements ³ (1) 24V 18Ah Lithium ion (standard) (1) 25.2V 10.4 Ah Lithium-lon (optional)	al)
Battery Weight ⁴ 7 lbs. (standard) 5.2 lbs. (optional)	
Battery Charger 2-amp, off-board (standard, optional)	
Color	

HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC OR DOES NOT MEET CRITERIA

Features

- A compact and lightweight design makes traveling a breeze
- Back seat pocket for extra storage while on-the-go
- Folds in just a few simple steps for effortless transport
- Includes convenient, under seat mesh storage bag, mesh cup holder and 40" lap belt
- Large foot platform to accommodate the user's comfort
- Specialty foam seat design for extreme comfort

Options

- Stealth® Simplicity Cushion 18"x18"
- Airline-approved travel battery

Accessories

XLR USB charger

Military Patches







Note: Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Jazzy® Passport and its lithium-ion battery are permitted.

- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%
- 2 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%
- 4 Battery weight may vary based on battery manufacturer.
- 5 Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power chair only. Overall measurements will vary based on seating and accessory selections.







