

Shower Stool

HCPCS: E0245



FIGURE A



FIGURE B

- Aluminum frame is lightweight, durable and corrosion proof
- Impact-resistant, composite seat
- Legs are height adjustable in 1" increments, and are crack proof and tarnish resistant
- Easy, tool-free assembly
- Made using recycled materials

PRODUCTS

Description

RTL12004KD	Bath Stool, KD, White,RTL,1 EA
------------	--------------------------------

SPECIFICATIONS

Description

Carton Shipping Weight	5 lbs.
Outside Legs	12.5" (W) x 12.5" (D)
Seat Dimensions	13" (W) x 13" (D)
Seat Height	13.5" - 21"
Warranty	Limited Lifetime
Weight	3 lbs.
Weight Capacity	300 lbs

ACCESSORIES

Description

12002T	Tips for 12002,4/pkg., BX
--------	---------------------------