

TranSport Aluminum Transport Chair

HCPCS: E1038



FIGURE A



FIGURE B

- Weighs 17.5 lbs without footrests
- Aluminum frame is lightweight and strong
- Frame magnets keep footrests secure against frame to prevent swinging open when the chair is carried
- Padded flip-back arms are removable
- Upholstery is comfortable, lightweight and easy-to-clean
- Two-way back release allows to push or pull lever
- The transport chair conveniently folds for travel and storage (Figure A)
- Back folds forward
- Comes with tool-free height adjustable swing-away footrests
- Aluminum rear wheel locks
- Weight Capacity: 250 lbs

PRODUCTS

Description

| | |
|------|--|
| TS19 | Tran Sport Chair,FlipBack-RemArms,1/cs |
|------|--|

SPECIFICATIONS

Description

| | |
|---|-------------------------------------|
| Back Height | 16" |
| Carton Shipping Weight | 19.5 lbs |
| Overall Height | 36" |
| Overall Length w/ Riggings | 38" |
| Overall Width | 24.5" (W) (open) 10" (W) (closed) |
| Seat Dimensions | 19" (W) x 16" (D) |
| Seat To Floor Height | 18" |
| Warranty | Limited Lifetime |
| Weight Capacity | 250 lbs |
| Weight Without Footrests With Footrests | 17.5 lbs 20 lbs |

PARTS

Description

| | |
|------------|------------------------------------|
| STDS871 | Hand Grip,Sentra,1/ea, EA |
| STDS1011 | Wheel Bearing, 8mm x 22mm 1/ea, EA |
| STDS1060 | Wheelchair Fork Bearing,1/ea, EA |
| STDS2A1800 | Rear Wheel, DFL, TS19,1/ea, EA |
| STDS2A1801 | REAR HUB CAP,DFL,1 PR, PR |
| STDS2A184L | BRAKE ASSY LEFT,DFL,1/EA, EA |
| STDS2A1811 | REAR AXLE BOLT,DFL,1 PR, PR |
| STDS2A1806 | FRONT AXLE BOLT,DFL,1 PR, PR |

| PARTS | Description |
|-------------|---|
| STDS2A1810 | FRONT WHEEL,DFL,TS19,1/ea, EA |
| TRB-1 | BackSupportReleaseLever./Trans,1/ea, EA |
| STDS1A2307 | Seat Upholstery,TS19,1/ea, EA |
| STDS2A184R | BRAKE ASSY RIGHT,DFL,1/EA, EA |
| STDS1A2308 | Back Upholstery,TS19,1/ea, EA |
| STDS001013 | Folding Cane Hindge TR18,TS19,1/ea, EA |
| STDS1A2305R | RightBackCaneAssy,TS19,1/ea, EA |
| STDS2A1803 | Cross Frame Assist (Blue) DFL,1/EA, EA |
| STDS2A1807 | SEAT RAIL PAD,DFL,NITRO,1/EA, EA |
| STDS1A2318R | Brake Assembly, Right, TS19, 1/ea, EA |
| STDS1A2309R | RightFootrestAssembly,TS19,1/ea, EA |
| STDS1069-4 | Hanger Pins, 4 per PKG, PKG |