

VivaLIFT![®]
POWER RECLINERS
By Pride Mobility Products[®] Corp.

THE MOST
COMFORTABLE
RECLINER
IN THE WORLD!

THE 4-MOTOR RECLINER

WITH
COUNTLESS
POSITION
COMBINATIONS



Shown in Badlands
Walnut Fabric





ATLAS PLUS



Wireless
Charging
Station*

When you need a moment for yourself, the **Atlas Plus Collection** is there for you. Customize your comfort while you watch the big game or read your favorite book. Enjoy a refreshing beverage and keep it close with the cup holder, and conveniently charge your phone with the wireless charger. Take advantage of those rare moments of solitude with the Atlas Plus Recliner.

* The wireless charger is only compatible with smart phones that have wireless charging capabilities. Additionally, no other smart devices, (i.e. smart watches) should be used with the wireless charger.





Pride Mobility[®] (US) 800-800-8586 | (CAN) 888-570-1113
pridemobility.com | prideliftchairs.com |    

ATLAS PLUS*







COLLECTION	ATLAS PLUS
MODEL	PLR-2985M
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift
SIZE	Medium
WEIGHT CAPACITY	375 lbs.
SEAT-TO-FLOOR	19.5"
SEAT DEPTH	21"
SEAT WIDTH	22"
TOP OF BACK TO SEAT	26"
FULLY PADDED CHAISE	Yes
FOOTREST EXTENSION	Standard
SEATING TYPE	Total Comfort
CHAIR HEIGHT	42"
DISTANCE FROM WALL	16"
OVERALL WIDTH	39"
CHAIR WEIGHT	164 lbs.
WARRANTY	7-Yr. Pro-rated

FUNCTIONS AND FEATURES:

-  POWER HEADREST
-  POWER LUMBAR
-  USB REMOTE
-  LITHIUM BATTERY BACKUP
-  FOOTREST EXTENSION
-  POWER BACKREST

POSITION GUIDE:

-  RELAXED
-  RECLINED
-  UPWARD LIFT
-  INFINITE LAY FLAT

STANDARD FABRIC: BADLANDS



MUSHROOM

WALNUT

STEEL

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

*Atlas Plus is available exclusively to Authorized Pride Brick & Mortar Providers. No internet sales or advertising of the Atlas Plus. Other conditions apply.



Pride Mobility® (US) 800-800-8586 | (CAN) 888-570-1113

pridemobility.com | prideliftchairs.com |    